

Description of the digital library files for TitaniumBase EV

File Name:

TitaniumBase EV DENTSPLY Implants2014-15_v1.0.dme

Version: #1

Date of release: 2015-11-16

Compatible with:

3Shape Dental Designer software version 2014 and 2015

Description of the digital library files

This digital library enables the 3Shape Dental System to create a sufficient output to produce the ceramic part onto the Titanium Base EV within the ASTRA TECH Implant System EV.

Restrictions:

- These library files are compatible with the 3Shape Dental System in the referred version above and may not be used with any other CAD/CAM system.
- This library including the TitaniumBase EV assortment is available for all markets accepting the CE-mark.

Products included in the library

	REF	Description
TitaniumBase EV	25924	TitaniumBase EV 3.0 – 1 mm
	25925	TitaniumBase EV 3.0 – 2 mm
	25926	TitaniumBase EV 3.6 – 1 mm
	25927	TitaniumBase EV 3.6 – 2 mm
	25928	TitaniumBase EV 4.2 – 1 mm
	25929	TitaniumBase EV 4.2 – 2 mm
	25930	TitaniumBase EV 4.8 – 1 mm
	25931	TitaniumBase EV 4.8 – 2 mm
	25932	TitaniumBase EV 5.4 – 1 mm
	25933	TitaniumBase EV 5.4 – 2 mm
ATLANTIS™ FLO	34593	OsseoSpeed™ EV 3.0 (P-01)
	34594	OsseoSpeed™ EV 3.6 (P-02)
	34595	OsseoSpeed™ EV 4.2 (P-03)
	34596	OsseoSpeed™ EV 4.8 (P-04)
	34597	OsseoSpeed™ EV 5.4 (P-05)
ATLANTIS™ Intra-Oral FLO	35243	OsseoSpeed™ EV 3.0 (IO-P-01)
	35244	OsseoSpeed™ EV 3.6 (IO-P-02)
	35245	OsseoSpeed™ EV 4.2 (IO-P-03)
	35246	OsseoSpeed™ EV 4.8 (IO-P-04)
	35247	OsseoSpeed™ EV 5.4 (IO-P-05)
ELOS Accurate® IO	IO 3B-A SA	ELOS Accurate® IO SA OsseoSpeed™ EV 3.0
	IO 3B-B SA	ELOS Accurate® IO SA OsseoSpeed™ EV 3.6
	IO 3B-C SA	ELOS Accurate® IO SA OsseoSpeed™ EV 4.2
	IO 3B-D SA	ELOS Accurate® IO SA OsseoSpeed™ EV 4.8
	IO 3B-E SA	ELOS Accurate® IO SA OsseoSpeed™ EV 5.4

All trademarks and company names are the property of their respective owners.