


## ATLANTIS™ ISUS clinic-to-laboratory prescription form

Clinician's name \_\_\_\_\_

Patient name \_\_\_\_\_



Implant position	Manufacturer/Company	Implant system	Implant level	Abutment system	Polished surface

### ATLANTIS™ ISUS Bar



- Titanium
- Cobalt-chrome

### ATLANTIS™ ISUS Bridge



- Titanium
- Cobalt-chrome

### ATLANTIS™ ISUS Hybrid



- Titanium
- Cobalt-chrome

### ATLANTIS™ ISUS 2in1



Only available in titanium

## Removable prostheses

---



### ATLANTIS™ ISUS 2in1

Provides primary and secondary suprastructures for a removable solution. The primary is fixed to implants, while the secondary attaches to the primary using friction, and additional retention elements.

The secondary can be a hybrid or bridge suprastructure finished with custom teeth and denture resin.

ATLANTIS™ ISUS 2in1 is only available in titanium.



### ATLANTIS™ ISUS Bar

Provides a suprastructures for removable dentures. Bars can be designed with extensions that follow the gingival profile. Moreover, they can include a combination of various attachment options. Both standard and custom-shaped bars are available, depending on case requirements.

## Fixed prostheses

---



### ATLANTIS™ ISUS Bridge

Provides a full-anatomical base for ceramic layering or composite layering techniques.



### ATLANTIS™ ISUS Hybrid

Provides a surface with retention elements that can be finished with custom teeth and denture resin.