

## ATLANTIS™ abutments implant compatibility chart



Manufacturer	Implant	ATLANTIS™ Abutment		ATLANTIS™ Crown		ATLANTIS™ Crown Abutment		ATLANTIS™ Conus Abutment	
		Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia	Titanium	Zirconia	Custom	Over-denture
								Titanium Gold-shaded ti	Titanium
<b>BioHorizons</b>	<b>INTERNAL</b>								
	Internal 3.5, 4.5, 5.7	x	x	x	x	x*	x	x	x*
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	x*	x	x	x*
	<b>EXTERNAL</b>								
	External implant 3.5	x**						x**	
External implant 4.0, 5.0	x**	x				x	x**		
<b>Biomet 3i</b>	<b>INTERNAL</b>								
	Certain™ 4/3, 5/4	x	x	x	x	x*	x	x	x*
	Certain™ 3.25, 4.0, 5.0, 6.0	x	x	x	x	x*	x	x	x*
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	x	x	x	x	x*	x	x	x*
	Certain™ XP 4/5, 5/6	x	x	x	x	x	x	x	x*
	<b>EXTERNAL</b>								
	Miniplant™ 3.25	x	x			x	x	x	x
Osseotite™ 3.75, 4.0, 5.0, 6.0	x	x			x	x	x	x	
Osseotite™ XP 4/5, 5/6	x	x			x	x	x	x	
<b>Camlog</b>	<b>INTERNAL</b>								
	Screw-line Implant 3.3	x		x		x		x	x
Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0	x	x	x	x	x	x	x	x	
<b>DENTSPLY Implants</b> <small>OsseoSpeed implants are part of the ASTRA TECH Implant System.</small>	<b>INTERNAL</b>								
	ANKYLOS® C/X – C/ 3.5, 4.5, 5.5, 7.0 (non-indexed)	x	x					x	
	ANKYLOS® C/X – /X 3.5, 4.5, 5.5, 7.0 (indexed)	x	x	x	x			x	x
	OsseoSpeed™ EV 3.0 – Green	x		x		x			
	OsseoSpeed™ EV 3.6 – Purple	x	x	x	x	x	x	x	x
	OsseoSpeed™ EV 4.2 – Yellow	x	x	x	x	x	x	x	x
	OsseoSpeed™ EV 4.8 – Blue	x	x	x	x	x	x	x	x
	OsseoSpeed™ EV 5.4 – Brown	x	x	x	x	x	x	x	x
	OsseoSpeed™ Profile EV 4.2 Yellow	x	x	x		x		x	x
	OsseoSpeed™ Profile EV 4.8 Blue	x	x	x		x		x	x
	OsseoSpeed™/OsseoSpeed™ TX 3.0 – Yellow	x		x		x			
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 – Aqua	x	x	x	x	x	x	x	x
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 – Lilac	x	x	x	x	x	x	x	x
	OsseoSpeed™ TX Profile 4.5, 5.0	x		x		x		x	x
	OsseoSpeed™ 4.0 – Blue (Previous)	x**						x**	
	XiVE® S 3.0	x		x		x			
	XiVE® S 3.4, 3.8, 4.5, 5.5	x	x	x		x	x	x	x
FRIALIT™ 3.4, 3.8, 4.5, 5.5	x	x	x		x	x	x	x	
<b>Keystone Dental</b>	<b>INTERNAL</b>								
	PrimaConnex™ SD 3.3, 3.5, Genesis™ 3.8	x	x	x	x	x	x	x	x
	PrimaConnex™ RD 4.0, 4.1, Genesis™ 4.5	x	x	x	x	x	x	x	x
PrimaConnex™ WD 5.0, Genesis™ 5.5, 6.5	x	x	x	x	x	x	x	x	

All trademarks and company names are the property of their respective owners.

\* 2 day extra Turnaround time.

\*\* 3 day extra Turnaround time.

Please turn over ➡

Manufacturer	Implant	ATLANTIS™ Abutment		ATLANTIS™ Crown		ATLANTIS™ Crown Abutment		ATLANTIS™ Conus Abutment	
		Titanium Gold-shaded ti	Zirconia	Titanium Gold-shaded ti	Zirconia	Titanium	Zirconia	Custom	Over-denture
								Titanium Gold-shaded ti	Titanium
Nobel Biocare	<b>INTERNAL</b>								
	NobelActive™ 3.0	x		x		x			
	NobelActive™ NP 3.5	x	x	x	x	x	x	x	x
	NobelActive™ RP 4.3, 5.0	x	x	x	x	x	x	x	x
	NobelReplace™ NP 3.5	x	x	x	x	x	x	x	x
	NobelReplace™ RP 4.3	x	x	x	x	x	x	x	x
	NobelReplace™ WP 5.0	x	x	x	x	x	x	x	x
	NobelReplace™ 6.0	x	x	x	x	x	x	x	x
	NobelReplace™ Conical Connection NP 3.5	x	x	x	x	x	x	x	x
	NobelReplace™ Conical Connection RP 4.3, 5.0	x	x	x	x	x	x	x	x
	<b>EXTERNAL</b>								
	Brånemark System™ NP 3.3	x				x		x	x
	Brånemark System™ RP 3.75, 4.0	x	x			x	x	x	x
Brånemark System™ WP 5.0, 5.5	x	x			x	x	x		
Straumann	<b>INTERNAL</b>								
	Bone Level 3.3 NC (Narrow CrossFit™)	x	x	x	x	x	x	x	x
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x	x	x	x	x	x	x	x
	Standard/Standard Plus 4.8 RN (Regular Neck)	x		x		x		x	x
	Standard/Standard Plus 4.8 WN (Wide Neck)	x		x		x		x	x
	<b>EXTERNAL</b>								
Standard Plus 3.3 NN (Narrow Neck)	x**						x**		
Zimmer Dental	<b>INTERNAL</b>								
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x	x			x	x	x	x
	Tapered Screw-Vent™ 3.5, 4.5, 5.7 friction fit	x				x		x	x
	Screw-Vent™ 3.3, 3.7, 4.5	x	x			x	x	x	x
	Screw-Vent™ 3.3, 3.7, 4.5 friction fit	x				x		x	x
SwissPlus™ 4.8	x				x		x	x	

All trademarks and company names are the property of their respective owners.

\* 2 day extra Turnaround time.

\*\* 3 day extra Turnaround time.